

Anti-6XHIS antibody (STJA0040603)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	WB
Host / Source	Chicken

PRODUCT PROPERTIES

Clonality	Polyclonal
Concentration	10 mg/mL
Conjugation	Unconjugated
Purification	Mixture of IgY fraction and affinity-purified antibodies
Dilution Range	WB 1:5000-1:10000
Formulation	Sodium phosphate-buffered (10 mM) isotonic (0.9%, w/v) saline (PBS, " pH 7.2) (50% by volume) with glycerol (50% by volume) as an anti-freezing agent and sodium azide (0.02%, w/v) as a preservative"
Isotype	IgY
Storage Instruction	Store at-20°C in the dark. Under these conditions, the antibodies should have a shelf life of at least twelve months, provided they remain sterile.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---