

Anti-Ki-67 antibody (STJA0040605)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ICC
Host / Source	Chicken
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot Dependent
Conjugation	Unconjugated
Purification	Total IgY fraction
Dilution Range	WB 1:5000 IHC 1:1000-1:5000 ICC 1:1000-1:5000
Formulation	PBS + 10 mM NaN3
Isotype	IgY
Molecular Weight	345-395 kDa
Storage Instruction	Recommended that the undiluted antibody be aliquoted into smaller working volumes (10-30 uL/vial depending on usage) upon arrival and stored long term at -20°C C or -80°C C, while keeping a working aliquot stored at 4°C C for short term. Avoid freeze/t

TARGET INFORMATION

Gene ID	4288
Gene Symbol	MKI67
UniProt ID	KI67_HUMAN
Immunogen Sequence	Human
Specificity	Specific for endogenous levels of the ~345 & 395 kDa Ki-67 isoforms. Additional bands are thought to be breakdown products of these 2 isoforms.

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



antibody WB analysis Ki-67 detection in HeLa cell lysate with anti-chicken HRP STJA0040605