

Anti-Neu-N antibody (STJA0040621)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ICC/IHC/WB
Host / Source	Chicken
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot dependent
Conjugation	Unconjugated
Purification	Total IgY fraction
Dilution Range	WB 1:250-1:500 WB Brain 1:250-1:500 IHC 1:2000-1:4000 ICC 1:2000-1:4000
Formulation	Concentrated IgY preparation plus 0.02% NaN3
Isotype	IgY
Molecular Weight	46, 48 kDa
Storage Instruction	Recommended that the undiluted antibody be aliquoted into smaller working volumes (10-30 uL/vial depending on usage) upon arrival and stored long term at -20°C C or -80°C C, while keeping a working aliquot stored at 4°C C for short term. Avoid freeze/t

TARGET INFORMATION


Gene ID	146713
Gene Symbol	RBFOX3
UniProt ID	RFOX3_HUMAN
Immunogen Sequence	Human
Specificity	Does not cross react with Fox1 or Fox2

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



IHC staining rat brain cerebellum using anti-antibody green neuronal labeling STJA0040621



IHC staining rat brain hippocampus with anti-antibody neuronal cells in green STJA0040621