

Anti-Perpherin antibody (STJA0040628)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC/ICC/WB
Host / Source	Chicken
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Concentration	2 mg/mL
Conjugation	Unconjugated
Purification	IgY Fraction
Dilution Range	WB 1:2000-1:5000 IHC 1:2000-1:5000 ICC 1:2000-1:5000
Formulation	Phosphate-buffered (10 mM) isotonic (0.9%, w/v) saline (PBS" pH 7.2) with bovine serum albumin (BSA, 0.5%) added to prevent absorption to the plastic and sodium azide (0.02%, w/v) added as a preservative"
Isotype	IgY
Molecular Weight	60 kDa
Storage Instruction	Store at 4°C in the dark. Under these conditions, the antibodies should have a shelf life of at least twelve months, provided they remain sterile. For longer term storage, aliquot and freeze to avoid freeze-thaw of the antibody.

TARGET INFORMATION

Gene ID	5630
Gene Symbol	PRPH
UniProt ID	PERI_HUMAN
Immunogen Sequence	Rat

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---