

Anti-Alpha Synuclein antibody (Full-Length) (STJA0040634)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC/ICC/WB
Host / Source	Chicken
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Concentration	10 mg/mL
Conjugation	Unconjugated
Purification	IgY Fraction
Dilution Range	WB 1:1000-1:5000 IHC 1:100-1:500 ICC 1:100-1:500
Formulation	Phosphate-buffered (10 mM) isotonic (0.9%, w/v) saline (PBS, " pH 7.2) with sodium azide (0.02%, w/v) added as a preservative."
Isotype	IgY
Molecular Weight	14 kDa
Storage Instruction	Store at 4°C in the dark. Under these conditions, the antibodies should have a shelf life of at least twelve months, provided they remain sterile. For longer term storage, aliquot and freeze to avoid freeze-thaw of the antibody.

TARGET INFORMATION

Gene ID	6622
Gene Symbol	SNCA
UniProt ID	SYUA_HUMAN
Immunogen Region	Full-Length
Immunogen Sequence	Human

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---