

Anti-Casein kinase I isoform alpha antibody (Full-Length) (STJA0040646)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ICC/IHC-F/WB
Host / Source	Chicken
Reactivity	Bovine/Chicken/Horse/Human/Mouse/Pig/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	IgY preparation
Dilution Range	WB 1:5000-1:10000 ICC 1:500-1:1000
Formulation	Lyophilized IgY preparation, with sodium azide.
Isotype	IgY
Molecular Weight	38 kDa
Storage Instruction	Spin vial briefly before opening. Reconstitute with 100 μ L sterile-filtered, ultrapure water. Centrifuge to remove any insoluble material. After reconstitution of lyophilized antibody, aliquot and store at -20°C for a higher stability. Avoid freeze-thaw.

TARGET INFORMATION

Gene ID	1452
Gene Symbol	CSNK1A1
UniProt ID	KC1A_HUMAN
Immunogen Region	Full-Length
Immunogen Sequence	Human
Specificity	Antibody recognizes CK1 by Western blot and immunohistochemistry. Cross reactivity.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---