

Anti-Microtubule-associated protein 2 antibody (STJA0040649)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ICC/IHC-F/IHC-P/WB
Host / Source	Chicken
Reactivity	Bovine/Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	IgY preparation
Dilution Range	WB 1:10000-1:20000 IHC 1:5000-1:10000 ICC 1:5000-1:10000
Formulation	Lyophilized IgY preparation, with sodium azide.
Isotype	IgY
Storage Instruction	Spin vial briefly before opening. Reconstitute with 50 μ L sterile-filtered, ultrapure water. Centrifuge to remove any insoluble material. After reconstitution of lyophilized antibody, aliquot and store at -20°C for a higher stability. Avoid freeze-thaw.

TARGET INFORMATION

Gene ID	4133
Gene Symbol	MAP2
UniProt ID	MTAP2_HUMAN
Immunogen Sequence	Bovine
Specificity	The specificity of this antibody has been confirmed by IC. Hu, rat, Ms, Bov.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---