

Anti-Neurofilament light polypeptide antibody (Full-Length) (STJA0040652)

GENERAL INFORMATION		
	Product Type	Primary antibodies
	Applications	IF/ICC/IHC/WB
	Host / Source	Chicken
	Reactivity	Bovine/Horse/Human/Mouse/Pig/Rat
PRODUCT PROPERTIES		
	Clonality	Polyclonal
	Conjugation	Unconjugated
	Purification	IgY preparation
	Dilution Range	WB 1:5000-1:10000 IHC 1:1000-1:5000 ICC 1:1000-1:5000
	Formulation	Lyophilized IgY preparation, with sodium azide.
	Isotype	IgY
	Storage Instruction	Spin vial briefly before opening. Reconstitute with 50 μL sterile-filtered, ultrapure water. Centrifuge to remove any insoluble material. Store lyophilized antibody at 2-8°C After reconstitution of lyophilized antibody, aliquot and store at-20°C for
TARGET INFORMATION		
	Gene ID	4747
	Gene Symbol	NEFL
	UniProt ID	NFL_HUMAN
	Immunogen Region	Full-Length
	Immunogen Sequence	Human
	Specificity	Species cross-reactivity includes human, rat, mouse, cow and pig. Predicted to react with other mammalian tissues due to sequence homology.
ADDITIONAL INFORMATION		
	Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.