

Anti-Peripherin antibody (Full-Length) (STJA0040655)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ICC/IHC-F/IHC-P/WB
Host / Source	Chicken
Reactivity	Cat/Human/Mouse/Other Mammals/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	IgY preparation
Dilution Range	WB 1:5000-1:10000 ICC 1:500-1:1000
Formulation	Lyophilized IgY preparation, with sodium azide.
Isotype	IgY
Molecular Weight	57 kDa
Storage Instruction	Spin vial briefly before opening. Reconstitute with 50 μ L sterile-filtered, ultrapure water. Centrifuge to remove any insoluble material. After reconstitution of lyophilized antibody, aliquot and store at -20°C for a higher stability. Avoid freeze-tha

TARGET INFORMATION

Gene Symbol	Prph Uniprot ID= PERL_RAT "
Immunogen Region	Full-Length
Immunogen Sequence	Rat
Specificity	The specificity of this antibody has been confirmed by WB. This antibody detects ~57 kDa Peripherin protein. A suitable control tissue is rat spinal cord or peripheral nerve homogenate. Hu, rat, Ms, Fel, and other mammals.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---