

Anti-C-reactive protein antibody (C-Term) (STJA0040660)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB
Host / Source	Chicken
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Affinity-purified
Dilution Range	WB 1:2000-1:5000
Formulation	PBS with 0.02% sodium azide.
Isotype	IgY
Storage Instruction	Short term storage at 2-8°C for one week. At -20°C as an undiluted liquid for up to 12 months.

TARGET INFORMATION

Gene ID	12944
Gene Symbol	Crp
UniProt ID	CRP_MOUSE
Immunogen Region	C-Term
Immunogen Sequence	Mouse
Specificity	C-reactive protein mouse.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---