

Anti-Leptin antibody (STJA0040663)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB
Host / Source	Chicken
Reactivity	Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Affinity-purified
Dilution Range	WB 1:1000-1:5000
Formulation	PBS with 0.02% sodium azide.
Isotype	IgY
Storage Instruction	Short term storage at 2-8°C for one week. At -20°C as an undiluted liquid for up to 12 months.

TARGET INFORMATION

Gene ID	25608
Gene Symbol	Lep
UniProt ID	LEP_RAT
Immunogen Sequence	Rat
Specificity	Rat Leptin.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---