

Anti-Brain-derived neurotrophic factor antibody (73-87 & 194-208 aa) (STJA0040665)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	ELISA
Host / Source	Chicken
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	IgY
Formulation	PBS with 0.02% Sodium Azide
Isotype	IgY
Storage Instruction	Short term storage at 2-8°C for one week. At -20°C as an undiluted liquid for up to 12 months.

TARGET INFORMATION

Gene ID	627
Gene Symbol	BDNF
UniProt ID	BDNF_HUMAN
Immunogen Region	73-87 & 194-208 aa
Immunogen Sequence	Human
Specificity	Mammalian BDNF.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---