

Anti-Glial cell line-derived neurotrophic factor antibody (101-118 & 199-211 aa) (STJA0040666)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	ELISA/ICC/IHC-F/WB
Host / Source	Chicken
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Affinity-purified
Dilution Range	WB 1:500-1:2000 IHC 1:50-1:100
Formulation	PBS with 0.02% sodium azide.
Isotype	IgY
Storage Instruction	Short term storage at 2-8°C for one week. At -20°C as an undiluted liquid for up to 12 months.

TARGET INFORMATION

Gene ID	2668
Gene Symbol	GDNF
UniProt ID	GDNF_HUMAN
Immunogen Region	101-118 & 199-211 aa
Immunogen Sequence	Human
Specificity	Human GDNF, cross reactivity to mouse, rat.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---