

**Anti-Appetite-regulating hormone, active antibody (24-33 aa) (STJA0040673)****GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	ELISA
Host / Source	Chicken
Reactivity	Human/Mouse/Rat

**PRODUCT PROPERTIES**

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Affinity-purified
Formulation	PBS with 0.02% sodium azide.
Isotype	IgY
Storage Instruction	Short term storage at 2-8°C for one week. At-20°C as an undiluted liquid for up to 12 months.

**TARGET INFORMATION**

Gene ID	<a href="#">51738</a>
Gene Symbol	<a href="#">GHRL</a>
UniProt ID	<a href="#">GHRL_HUMAN</a>
Immunogen Region	24-33 aa
Immunogen Sequence	Human
Specificity	Active Ghrelin (24-33) : Affinity purified, human, mouse specific.

**ADDITIONAL INFORMATION**

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---