

Anti-Ciliary neurotrophic factor antibody (11-25 & 186-200 aa) (STJA0040674)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	ELISA/WB
Host / Source	Chicken
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Affinity-purified
Dilution Range	WB 1:500-1:2000 IHC 1:50-1:500
Formulation	PBS with 0.02% sodium azide.
Isotype	IgY
Storage Instruction	Short term storage at 2-8°C for one week. At -20°C as an undiluted liquid for up to 12 months.

TARGET INFORMATION

Gene ID	1270
Gene Symbol	CNTF
UniProt ID	CNTF_HUMAN
Immunogen Region	11-25 & 186-200 aa
Immunogen Sequence	Human
Specificity	CNTF, Cross reactivity to other species is expected.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---