

Anti-Fatty acid-binding protein, adipocyte antibody (71-82 aa+97-111 aa) (STJA0040676)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	ELISA/WB
Host / Source	Chicken
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Affinity-purified
Dilution Range	WB 1:1000-1:5000
Formulation	PBS with 0.02% sodium azide.
Isotype	IgY
Molecular Weight	16 kDa
Storage Instruction	Short term storage at 2-8°C for one week. At-20°C as an undiluted liquid for up to 12 months.

TARGET INFORMATION

Gene ID	2167
Gene Symbol	FABP4
UniProt ID	FABP4_HUMAN
Immunogen Region	71-82 aa+97-111 aa
Immunogen Sequence	Human
Specificity	This antibody detects Fatty acid-binding protein, adipocyte in various mouse and rat fat cell lysates at approx 16 kDa. Human, mouse, rat.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---