

Anti-Orexin B antibody (N-Term) (STJA0040680)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA
Host / Source	Chicken
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Affinity-purified
Formulation	PBS with 0.02% sodium azide.
Isotype	IgY
Storage Instruction	Short term storage at 2-8°C for one week. At -20°C as an undiluted liquid for up to 12 months.

TARGET INFORMATION

Gene ID	3060
Gene Symbol	HCRT
UniProt ID	OREX_HUMAN
Immunogen Region	N-Term
Immunogen Sequence	Human
Specificity	Human Orexin B (1-14) : Affinity purified.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---