

**Anti-Peroxisome proliferator-activated receptor gamma antibody (115-129 aa) (STJA0040681)****GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	WB
Host / Source	Chicken
Reactivity	Human

**PRODUCT PROPERTIES**

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Affinity-purified
Formulation	PBS with 0.02% sodium azide.
Isotype	IgY
Storage Instruction	Short term storage at 2-8°C for one week. At-20°C as an undiluted liquid for up to 12 months.

**TARGET INFORMATION**

Gene ID	<a href="#">5468</a>
Gene Symbol	<a href="#">PPARG</a>
UniProt ID	<a href="#">PPARG_HUMAN</a>
Immunogen Region	115-129 aa
Immunogen Sequence	Human
Specificity	Human Peroxisome proliferator-activated receptor gamma.

**ADDITIONAL INFORMATION**

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---