

Anti-Mitochondrial uncoupling protein 3 antibody (C-Term) (STJA0040690)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB
Host / Source	Chicken
Reactivity	Human/Mouse/Pig/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Affinity-purified
Dilution Range	WB 1:1000-1:5000
Formulation	PBS with 0.02% sodium azide.
Isotype	IgY
Molecular Weight	35 kDa
Storage Instruction	Short term storage at 2-8°C for one week. At -20°C as an undiluted liquid for up to 12 months.

TARGET INFORMATION

Gene ID	7352
Gene Symbol	UCP3
UniProt ID	UCP3_HUMAN
Immunogen Region	C-Term
Immunogen Sequence	Human
Specificity	This antibody detects UCP3 at approx 35 kDa. Human, mouse, rat, Porcine.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---