

Anti-Growth associated protein 43 antibody (C-Term) (STJA0040693)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ICC/WB
Host / Source	Chicken
Reactivity	Bovine/Human/Other Mammals/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Concentrated purified IgY
Dilution Range	WB 1:1000-1:5000 ICC 1:500-1:1000
Formulation	Lyophilized IgY preparation, with sodium azide.
Isotype	IgY
Molecular Weight	43 kDa
Storage Instruction	Spin vial briefly before opening. Reconstitute with 100 μ L sterile-filtered, ultrapure water. Centrifuge to remove any insoluble material. Store lyophilized, unopened vial at 2-8°C or lower. After reconstitution, prepare aliquots and store at -20°C to

TARGET INFORMATION

Gene ID	29423
Gene Symbol	Gap43
UniProt ID	NEUM_RAT
Immunogen Region	C-Term
Immunogen Sequence	Rat
Specificity	The antibody reacts with a 43 kDa band by Western blot on bovine hippocampus lysate. It has also been used successfully for immunocytochemistry. Human, rat and bovine. It is expected that it will work on other mammal tissues due to amino acid sequenc

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---