

Anti-GDNF family receptor alpha-1 antibody (STJA0040695)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ICC/WB
Host / Source	Chicken
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Affinity-purified
Dilution Range	WB 1:5000-1:10000 ICC 1:500-1:2000
Formulation	Concentrated ammonium sulphate in PBS pH 7.4 containing no preservatives.
Isotype	IgY
Storage Instruction	Store antibody stock solution at +4°C upon receipt. DO NOT FREEZE. This product is an (NH ₄) ₂ SO ₄ precipitate. To use, mix well by pipetting or vortexing gently prior use, then withdraw the desired amount of antibody stock solution and dilute into you

TARGET INFORMATION

Gene ID	2674
Gene Symbol	GFRA1
UniProt ID	GFRA1_HUMAN
Immunogen Sequence	Human
Specificity	Human GFR alpha-1. No cross reaction with hGFR alpha-2, hGFR alpha-3, hGFR alpha-4.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---