

Anti-Microtubule-associated protein tau antibody (Full-Length) (STJA0040700)

GENERAL INFORMATION	
Product Type	Primary antibodies
Applications	ICC/WB
Host / Source	Chicken
Reactivity	Bovine/Chicken/Horse/Human/Mouse/Pig/Rat
PRODUCT PROPERTIES	
Clonality	Polyclonal
Conjugation	Unconjugated
Purification	IgY preparation
Dilution Range	WB 1:5000-1:10000 ICC 1:500-1:1000
Formulation	Lyophilized IgY preparation, with sodium azide.
Isotype	IgY
Molecular Weight	48 kDa to 67 kDa, with an additional band at 100 kDa
Storage Instruction	Spin vial briefly before opening. Reconstitute with 100 μ L sterile-filtered, ultrapure water. Centrifuge to remove any insoluble material. Store lyophilized, unopened vial at 2-8°C or lower. After reconstitution, prepare aliquots and store at -20°C to
TARGET INFORMATION	
Gene ID	4137
Gene Symbol	MAPT
UniProt ID	TAU_HUMAN
Immunogen Region	Full-Length
Immunogen Sequence	Human
Specificity	The antibody reacts with multiple closely spaced bands covering the region of the blot from 48 kDa to 67 kDa, with an additional band at 100 kDa. It has also been used successfully for immunocytochemistry. Expected to react with horse, cow, pig, chic
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.