

Anti-Casein, A1 beta antibody (78-86) (STJA0040704)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	ELISA/WB
Host / Source	Chicken
Reactivity	Bovine

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Affinity-purified
Dilution Range	WB 1:1000-1:5000
Formulation	PBS, pH 7.4, containing 0.02% sodium azide as preservative.
Isotype	IgY
Storage Instruction	Maintain unopened vial at -20°C for up to 12 months after date of receipt. After opening maintain at -20°C in undiluted aliquots for up to 6 months. For short-term storage, keep aliquot at 2-8°C for up to one week. Avoid repeated freeze-thaw cycles.

TARGET INFORMATION

Immunogen Region	78-86
Immunogen Sequence	Bovine
Specificity	The antibody is specific to A1 beta casein by Western blot. No cross-reactivity with A2 beta casein is seen. Species cross-reactivity not tested.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---