

Anti-Alpha synuclein antibody (61-95) (STJA0040705)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ICC/IHC-F/WB
Host / Source	Chicken
Reactivity	Bovine/Chicken/Horse/Human/Mouse/Pig/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Affinity-purified
Dilution Range	WB 1:1000-1:2000 IHC 1:1000-1:2000 ICC 1:1000-1:2000
Formulation	PBS, pH 7.4, containing 0.02% sodium azide as preservative.
Isotype	IgY
Storage Instruction	Spin vial briefly before opening. Reconstitute with 50 μ L sterile-filtered, ultrapure water to achieve a 1 mg/mL concentration. Centrifuge to remove any insoluble material. Store lyophilized antibody at 2-8°C. After reconstitution divide into aliquot

TARGET INFORMATION

Gene ID	6622
Gene Symbol	SNCA
UniProt ID	SYUA_HUMAN
Immunogen Region	61-95
Immunogen Sequence	Human
Specificity	Reacts with human, horse, cow, pig, chicken, rat, mouse.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---