

Anti-Calbindin antibody (Full-Length) (STJA0040706)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-F/WB
Host / Source	Chicken
Reactivity	Bovine/Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	IgY preparation
Dilution Range	WB 1:1000-1:5000 IHC 1:1000-1:5000
Formulation	Lyophilized IgY preparation, with sodium azide.
Isotype	IgY
Storage Instruction	Spin vial briefly before opening. Reconstitute with 50 μ L sterile-filtered, ultrapure water. Centrifuge to remove any insoluble material. Store lyophilized antibody at 2-8°C. After reconstitution divide into aliquots and store at 20°C for long-term s

TARGET INFORMATION

Gene ID	793
Gene Symbol	CALB1
UniProt ID	CALB1_HUMAN
Immunogen Region	Full-Length
Immunogen Sequence	Human
Specificity	Reacts with human, cow, rat, mouse. Antibody is specific for calbindin and does not recognize closely related proteins parvalbumin and calretinin as determined by Western Blotting.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---