

## Anti-Ki-67 antibody (1-300 & 1, 111-1, 490) (STJA0040709)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ICC/WB
Host / Source	Chicken
Reactivity	Human

### PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	IgY preparation
Dilution Range	WB 1:2000-1:5000 ICC 1:1000-1:5000
Formulation	Lyophilized IgY preparation, with sodium azide.
Isotype	IgY
Storage Instruction	Spin vial briefly before opening. Reconstitute with 50 $\mu$ L sterile-filtered, ultrapure water. Centrifuge to remove any insoluble material. Store lyophilized antibody at 2-8°C. After reconstitution divide into aliquots and store at 20°C for long-term s

### TARGET INFORMATION

Gene ID	<a href="#">4288</a>
Gene Symbol	<a href="#">MKI67</a>
UniProt ID	<a href="#">K167_HUMAN</a>
Immunogen Region	1-300 & 1, 111-1, 490
Immunogen Sequence	Human
Specificity	Reacts with human only. Does not react with mouse or rat.

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---