

Anti-Apolipoprotein E antibody (STJA0040717)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA
Host / Source	Goat
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Affinity-purified
Formulation	Lyophilized from a solution containing 50 mM Tris, pH 7.5, 0.4 M NaCl, 0.01 M EDTA, 0.1% trehalose, 0.07% sodium azide.
Isotype	IgG
Storage Instruction	Spin vial briefly before opening. Reconstitute in 100 μ L sterile-filtered, ultrapure water. Centrifuge to remove any insoluble material. Store lyophilized antibody at 2-8°C. After reconstitution keep aliquots at -20°C to -80°C for higher stability. Avo

TARGET INFORMATION

Gene ID	348
Gene Symbol	APOE
UniProt ID	APOE_HUMAN
Immunogen Sequence	Human
Specificity	Human. Species cross-reactivity not tested.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---