

Anti-Green fluorescent protein antibody (STJA0040719)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ICC/WB
Host / Source	Goat
Reactivity	Species Independent

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Affinity-purified
Dilution Range	WB 1:1000-1:5000 ICC 1:2000-1:5000
Formulation	Lyophilized from PBS buffer pH 7.2-7.6 with 0.1% trehalose, and sodium azide
Isotype	IgG
Molecular Weight	27 kDa
Storage Instruction	Spin vial briefly before opening. Reconstitute with 100 μ L sterile-filtered, ultrapure water to achieve a 1 mg/mL concentration. Centrifuge to remove any insoluble material. Store lyophilized antibody at 2-8°C. After reconstitution divide into aliquot

TARGET INFORMATION

Immunogen Sequence	Jellyfish
Specificity	Antibody recognizes GFP protein in GFP-transfected HEK293 cells, but not in non-transfected control cells. Specific for GFP, does not cross-react with mCherry.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---