

Anti-Melanocyte-stimulating hormone, alpha antibody (POMC) (STJA0040721)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-F
Host / Source	Guinea Pig
Reactivity	Human/Rat/Sheep

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Whole serum
Dilution Range	IHC 1.4305555555555556
Isotype	Mixed
Storage Instruction	Spin vial briefly before opening. Reconstitute in 50 μ L sterile-filtered, ultrapure water. Centrifuge to remove any insoluble material. Store lyophilized antibody at 2-8°C or lower. After reconstitution, divide into aliquots and store at -20°C, or low

TARGET INFORMATION

Gene ID	5443
Gene Symbol	POMC
UniProt ID	COLI_HUMAN
Immunogen Region	POMC
Immunogen Sequence	Human
Specificity	Specificity was demonstrated by immunohistochemistry. This antibody is known to react with human, rat and sheep. Other species have not yet been tested.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---