

## Anti-Neuropeptide RFRP-3 antibody (RFRP-1) (STJA0040726)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-F
Host / Source	Guinea Pig
Reactivity	Human/Rat/Sheep

### PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Whole serum
Dilution Range	IHC 1.4305555555555556
Isotype	Mixed
Storage Instruction	Spin vial briefly before opening. Reconstitute in 50 $\mu$ L sterile-filtered, ultrapure water. Centrifuge to remove any insoluble material. Store lyophilized antibody at 2-8°C or lower. After reconstitution, divide into aliquots and store at -20°C, or low

### TARGET INFORMATION

Gene ID	<a href="#">64111</a>
Gene Symbol	<a href="#">NPVF</a>
UniProt ID	<a href="#">NPVF_HUMAN</a>
Immunogen Region	RFRP-1
Immunogen Sequence	Human
Specificity	Specificity was demonstrated by immunohistochemistry. This antibody is known to react with rat and sheep.

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---