

**Anti-NGFR/p75 neurotrophin receptor antibody [ME20.4] (STJA0040730)****GENERAL INFORMATION**

|               |                                |
|---------------|--------------------------------|
| Product Type  | Primary antibodies             |
| Applications  | IHC                            |
| Host / Source | Mouse                          |
| Reactivity    | Cat/Dog/Human/Pig/Rabbit/Sheep |

**PRODUCT PROPERTIES**

|                     |                                                                                                                                                                                                                                                              |
|---------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Clonality           | Monoclonal                                                                                                                                                                                                                                                   |
| Conjugation         | Unconjugated                                                                                                                                                                                                                                                 |
| Purification        | IgG                                                                                                                                                                                                                                                          |
| Dilution Range      | IHC 2 µg/mL                                                                                                                                                                                                                                                  |
| Isotype             | IgG1                                                                                                                                                                                                                                                         |
| Storage Instruction | Spin vial briefly before opening. Reconstitute with the addition of 100 µL sterile-filtered, ultrapure water. Centrifuge to remove any insoluble material. After reconstitution keep aliquots at -20°C for a higher stability, and at 2-8°C with an appropri |

**TARGET INFORMATION**

|                    |                                                                                                                                                            |
|--------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Gene ID            | <a href="#">4804</a>                                                                                                                                       |
| Gene Symbol        | <a href="#">NGFR</a>                                                                                                                                       |
| UniProt ID         | <a href="#">TNR16_HUMAN</a>                                                                                                                                |
| Immunogen Sequence | Human                                                                                                                                                      |
| Specificity        | This antibody recognizes p75NTR (low affinity neurotrophin receptor) Reacts with human, cat, dog, pig, rabbit and sheep. Does not react with rat or mouse. |

**ADDITIONAL INFORMATION**

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|