

**Anti-NGFR/p75 neurotrophin receptor antibody [ME20.4] {FITC Ex: 490 nm, Em: 525 nm}
(STJA0040731)**

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC/IHC-F
Host / Source	Mouse
Reactivity	Cat/Dog/Human/Pig/Rabbit/Sheep

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	FITC Ex: 490 nm, Em: 525 nm
Purification	IgG
Dilution Range	IHC 2 µg/mL
Formulation	PBS pH 7.2-7.6
Isotype	IgG1
Storage Instruction	The antibody conjugate can be stored at 2-8°C for up to 4 months with the addition of appropriate antibacterial agent.

TARGET INFORMATION

Gene ID	4804
Gene Symbol	NGFR
UniProt ID	TNFR16_HUMAN
Immunogen Sequence	Human
Specificity	This antibody recognizes p75NTR (low affinity neurotrophin receptor) Reacts with human, cat, dog, pig, rabbit and sheep. Does not react with rat or mouse.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---