

**Anti-NGFR/p75 neurotrophin receptor antibody [MC192] {ATTO 488 Ex: 500 nm, Em: 520 nm}
(STJA0040732)****GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	FC/ICC/IHC-F/IHC-P
Host / Source	Mouse
Reactivity	Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	ATTO 488 Ex: 500 nm, Em: 520 nm
Purification	IgG
Dilution Range	IHC 2-5 µg/mL ICC 2-5 µg/mL
Isotype	IgG1
Storage Instruction	The antibody conjugate can be stored at 2-8°C for up to 4 months with the addition of appropriate antibacterial agent.

TARGET INFORMATION

Gene ID	24596
Gene Symbol	Ngfr
UniProt ID	TNR16_RAT
Immunogen Sequence	Rat
Specificity	MC192 recognizes the extracellular domain of the neurotrophin receptor p75NTR in rat and does not react with human or mouse NGFR. Reacts with rat. Does not react with mouse or human NGFR.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---