

Anti-Influenza hemagglutinin antibody (YPYDVPDYA) [HA.C5] (STJA0040738)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	ICC/IHC-F/WB
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity-purified
Dilution Range	WB 0.7361111111111116 IHC 0.1805555555555555 ICC 0.1805555555555555 IP 0.1805555555555555
Formulation	PBS, pH 7.4 with 0.05% sodium azide
Isotype	IgG
Storage Instruction	Stable for 1 year at -20°C from the date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to opening the cap. Aliquot to avoid repeated freezing and thawing.

TARGET INFORMATION

Immunogen Region	YPYDVPDYA
Immunogen Sequence	Human
Specificity	Antibody detects HA-tagged proteins (on amino- and carboxy termini) in transfected mammalian cells.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---