

Anti-Myc proto-oncogene protein antibody (EQKLISEEDL) [Myc.A7] (STJA0040740)

GENERAL INFORMATION	
Product Type	Primary antibodies
Applications	ICC/WB
Host / Source	Mouse
Reactivity	Species Independent
PRODUCT PROPERTIES	
Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity-purified
Dilution Range	WB 0.7361111111111116 IHC 0.1805555555555555 ICC 0.1805555555555555 IP 0.1805555555555555
Formulation	PBS, pH 7.4 with 0.05% sodium azide
Isotype	IgG1
Storage Instruction	Stable for 1 year at -20°C from the date of receipt. For maximum recovery of product, centrifuge the original vial after thawing and prior to opening the cap. Aliquot to avoid repeated freezing and thawing.
TARGET INFORMATION	
Gene ID	4609
Gene Symbol	MYC
UniProt ID	MYC_HUMAN
Immunogen Region	EQKLISEEDL
Immunogen Sequence	Human
Specificity	Antibody detects over-expressed proteins containing Myc epitope tag fused to either amino-or carboxy-termini of targeted proteins in transfected mammalian cells.
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.