

Anti-Simian Virus Type 5 tag antibody (GKIPNPLLGLDST) [V5.E10] (STJA0040742)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	ELISA/IHC-F/WB
Host / Source	Mouse
Reactivity	Species Independent

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity-purified
Dilution Range	WB 1:500-1:1000 IHC 1:50-1:100
Formulation	PBS, pH 7.4 with 0.05% sodium azide
Isotype	IgG
Storage Instruction	Stable for 1 year at -20°C from the date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to opening the cap. Aliquot to avoid repeated freezing and thawing.

TARGET INFORMATION

Immunogen Region	GKIPNPLLGLDST
Immunogen Sequence	Virus
Specificity	Antibody detects V5-tagged fusion protein.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---