

Anti-rRNA 2 O-methyltransferase fibrillarin antibody [38F3] (STJA0040745)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ICC/WB
Host / Source	Mouse
Reactivity	C.elegans/Drosophila/Human/Other Mammals/Rat/Yeast

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	IgG
Dilution Range	WB 1:500-1:10000 (see notes) ICC 1:500-1:5000 (see notes)
Formulation	Lyophilized hybridoma cell culture media with sodium azide
Isotype	IgG
Molecular Weight	34 kDa
Storage Instruction	Spin the vial briefly before opening it. Reconstitute in 250 uL sterile-filtered, room temperature, ultrapure water. Let the vial rehydrate for 5 minutes. Mix gently with a pipette tip to help the redissolution of the material. Centrifuge to remove a

TARGET INFORMATION

Immunogen Sequence	Yeast
Specificity	The specificity of this antibody has been confirmed by WB. This clone was selected because it is specific for the ~34 kDa Fibrillarin. It stains a single band on Western blotting and shows a clear and strong punctate staining of yeast nuclei. It can

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---