

Anti-Internexin alpha antibody [1D2] (STJA0040748)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ICC/IHC-F/IHC-P/WB
Host / Source	Mouse
Reactivity	Cat/Human/Mouse/Other Mammals/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	IgG
Dilution Range	WB 1:5000-1:10000 IHC 1:1000-1:5000 ICC 1:250-1:500
Formulation	Lyophilized from PBS buffer pH 7.2-7.6 with 0.1% trehalose, and sodium azide
Isotype	IgG1
Molecular Weight	64-66 kDa
Storage Instruction	Spin vial briefly before opening. Reconstitute with 250 μ L sterile-filtered, ultrapure water to achieve a 1 mg/mL concentration. Centrifuge to remove any insoluble material. After reconstitution of lyophilized antibody, aliquot and store at -20°C for

TARGET INFORMATION

Gene Symbol	Ina Uniprot ID= AINX_RAT "
Immunogen Sequence	Rat
Specificity	The specificity of this antibody has been confirmed by WB. This antibody is specific for the 64-66 kDa alpha-internexin protein. Molecular weight will depend on species. Hu, rat, Ms, Fel, and other mammals.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---