

Anti-Peripherin antibody (Full-Length) [7C5] (STJA0040754)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ICC/IHC-F/WB
Host / Source	Mouse
Reactivity	Cat/Human/Mouse/Other Mammals/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity-purified
Dilution Range	WB 0.5-2 μ g/mL IHC 1-5 μ g/mL ICC 1-5 μ g/mL
Formulation	Lyophilized from PBS buffer pH 7.2-7.6 with 0.1% trehalose, and sodium azide
Isotype	IgG
Molecular Weight	57 kDa
Storage Instruction	Spin vial briefly before opening. Reconstitute with 100 μ L sterile-filtered, ultrapure water to achieve a 1 mg/mL concentration. Centrifuge to remove any insoluble material. Store lyophilized, unopened vial at 2-8°C or lower. After reconstitution, pr

TARGET INFORMATION

Gene Symbol	Prph Uniprot ID= PERL_RAT "
Immunogen Region	Full-Length
Immunogen Sequence	Rat
Specificity	The specificity of this antibody has been confirmed by WB. This antibody detects ~57 kDa Peripherin protein. Human, mouse, feline. Predicted to react with other mammalian tissue.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---