

Anti-Polyubiquitin-B antibody [Ubi-1] (STJA0040757)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-F/IHC-P/WB
Host / Source	Mouse
Reactivity	Bovine/C.elegans/Chicken/Drosophila/Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	IgG
Dilution Range	WB 1:500-1:1000 IHC 0.388888888888889
Formulation	Lyophilized from PBS buffer pH 7.2-7.6 with 0.1% trehalose, and sodium azide
Isotype	IgG1
Molecular Weight	8.5 kDa
Storage Instruction	Spin vial briefly before opening. Reconstitute with 100 μ L sterile-filtered, ultrapure water to achieve a 1 mg/mL concentration. Centrifuge to remove any insoluble material. After reconstitution of lyophilized antibody, aliquot and store at -20°C for

TARGET INFORMATION

Gene ID	7314
Gene Symbol	UBB
UniProt ID	UBB_HUMAN
Immunogen Sequence	Human
Specificity	The specificity of this antibody has been confirmed by WB. This antibody detects ~8.5 kDa Ubiquitin. Hu, Bov, Chk, Drosophila, and C. elegans.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---