

Anti-Ubiquitin carboxyl-terminal hydrolase isozyme L1 antibody (Full-Length) [BH7] (STJA0040758)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ICC/IHC-F/WB
Host / Source	Mouse
Reactivity	Bovine/Human/Mouse/Pig/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	IgG
Dilution Range	WB 1:10000-1:20000 IHC 1:1000-1:5000 ICC 1:1000-1:5000
Formulation	Lyophilized from PBS buffer pH 7.2-7.6 with 0.1% trehalose, and sodium azide
Isotype	IgG1
Molecular Weight	24 kDa
Storage Instruction	Spin vial briefly before opening. Reconstitute with 100 μ L sterile-filtered, ultrapure water to achieve a 1 mg/mL concentration. Centrifuge to remove any insoluble material. After reconstitution of lyophilized antibody, aliquot and store at -20°C for

TARGET INFORMATION

Gene ID	7345
Gene Symbol	UCHL1
UniProt ID	UCHL1_HUMAN
Immunogen Region	Full-Length
Immunogen Sequence	Human
Specificity	The specificity of this antibody has been confirmed by WB. This antibody detects ~24 kDa UCHL1 enzyme. Suitable control tissue is rat spinal cord, brain, SHSY-5Y or HEK293 cell extract. Hu, rat, Bov.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---