

Anti-Lipid phosphate phosphohydrolase 3 antibody (179-196 aa) [7H7D3] (STJA0040759)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	FC/WB
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity-purified
Dilution Range	WB 2-3 µg/mL IHC 0.5-2 µg/mL
Formulation	Lyophilized from PBS pH 7.4
Isotype	Mix of IgG1, IgG2a & IgG2b
Storage Instruction	Spin vial briefly before opening. Reconstitute in 100 uL sterile-filtered, ultrapure water. Centrifuge to remove any insoluble material. At least 12 months after purchase at 2-8°C (lyophilized formulations). After reconstitution, aliquot and store at

TARGET INFORMATION

Gene ID	8613
Gene Symbol	PLPP3
UniProt ID	PLPP3_HUMAN
Immunogen Region	179-196 aa
Immunogen Sequence	Human
Specificity	Confirmed by over-expression of human LPP3 cDNA. Human.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---