

Anti-14-3-3 protein eta antibody (Full-Length) [3G12] (STJA0040760)

GENERAL INFORMATION	
Product Type	Primary antibodies
Applications	ICC/IHC-F/WB
Host / Source	Mouse
Reactivity	Human/Mouse/Rat
PRODUCT PROPERTIES	
Clonality	Monoclonal
Conjugation	Unconjugated
Purification	IgG
Dilution Range	WB 1:1000-1:5000 IHC 1:500-1:1000 ICC 1:500-1:1000
Formulation	Lyophilized from PBS buffer pH 7.2-7.6 with 0.1% trehalose, and sodium azide
Isotype	IgG
Storage Instruction	Spin vial briefly before opening. Reconstitute with 100 μ L sterile-filtered, ultrapure water to achieve a 1 mg/mL concentration. Centrifuge to remove any insoluble material. After reconstitution of lyophilized antibody, aliquot and store at -20°C for
TARGET INFORMATION	
Gene ID	7533
Gene Symbol	YWHAH
UniProt ID	1433F_HUMAN
Immunogen Region	Full-Length
Immunogen Sequence	Human
Specificity	Human 14-3-3 ETA protein.
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.