

Anti-Galectin-3 antibody (Full-Length) [5C21] (STJA0040762)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ICC/WB
Host / Source	Mouse
Reactivity	Human/Mouse/Other Mammals/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity-purified
Dilution Range	WB 1.430555555555556 ICC 0.7361111111111116
Formulation	Lyophilized from PBS buffer pH 7.2-7.6 with 0.1% trehalose, and sodium azide
Isotype	IgG1
Storage Instruction	Spin vial briefly before opening. Reconstitute with 100 μ L sterile-filtered, ultrapure water to achieve a 1 mg/mL concentration. Centrifuge to remove any insoluble material. After reconstitution of lyophilized antibody, aliquot and store at -20°C for

TARGET INFORMATION

Gene ID	3958
Gene Symbol	LGALS3
UniProt ID	LEG3_HUMAN
Immunogen Region	Full-Length
Immunogen Sequence	Human
Specificity	Human; mouse; rat; predicted to react with other species due to sequence homology;.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---