

Anti-Microtubule Associated Protein 2 antibody [5H11] (STJA0040766)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	ICC/IHC-F/WB
Host / Source	Mouse
Reactivity	Bovine/Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	IgG
Dilution Range	WB 1:5000-1:10000 IHC 1:1000-1:5000 ICC 1:1000-1:5000
Formulation	Lyophilized from PBS buffer pH 7.2-7.6 with 0.1% trehalose, and sodium azide
Isotype	IgG
Storage Instruction	Spin vial briefly before opening. Reconstitute with 100 μ L sterile-filtered, ultrapure water to achieve a 1 mg/mL concentration. Centrifuge to remove any insoluble material. At least 12 months after purchase at 2-8°C (lyophilized formulations). After

TARGET INFORMATION

Immunogen Sequence	Bovine
Specificity	The specificity of this antibody has been confirmed by WB and IHC against the antigen. Human; rat; mouse;.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---