

Anti-Visinin-like protein 1 antibody [2D11] (STJA0040773)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ICC/IHC-F/WB
Host / Source	Mouse
Reactivity	Bovine/Human/Mouse/Other Mammals/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	IgG
Dilution Range	WB 1:500-1:1000 IHC 1:500-1:1000 ICC 1:500-1:1000
Formulation	Lyophilized from PBS buffer pH 7.2-7.6 with 0.1% trehalose, and sodium azide
Isotype	IgG1
Molecular Weight	22 kDa
Storage Instruction	Spin vial briefly before opening. Reconstitute with 100 μ L sterile-filtered, ultrapure water to achieve a 1 mg/mL concentration. Centrifuge to remove any insoluble material. After reconstitution of lyophilized antibody, aliquot and store at -20°C for

TARGET INFORMATION

Gene ID	7447
Gene Symbol	VSNL1
UniProt ID	VISL1_HUMAN
Immunogen Sequence	Human
Specificity	The antibody reacts with a 22 kDa band by Western blot on bovine cerebellum lysate. It has also been used successfully for immunocytochemistry.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---