

Anti-Neurturin antibody [1B11] (STJA0040774)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	IgG
Dilution Range	WB 1:5000-1:10000
Formulation	PBS pH 7.4, with 0.1% sodium azide
Isotype	IgG1
Storage Instruction	Store at -20°C or -70°C upon receipt. After opening divide antibody into smaller aliquots and store at -20°C or -70°C for up to six months. Avoid multiple freeze-thaw cycles as product degradation may result.

TARGET INFORMATION

Gene ID	4902
Gene Symbol	NRTN
UniProt ID	NRTN_HUMAN
Immunogen Sequence	Human
Specificity	Human NTN (Neurturin).

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---