

Anti-Mesencephalic astrocyte-derived neurotrophic factor antibody (STJA0040775)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	IgG
Dilution Range	WB 0.5-2 µg/mL IHC 1-10 µg/mL ICC 1-10 µg/mL
Formulation	PBS pH 7.4, with 0.1% sodium azide
Isotype	IgG1
Storage Instruction	Store at -20°C or -70°C upon receipt. After opening maintain undiluted in smaller aliquots at -20°C or -70°C for up to 6 months. Avoid multiple freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7873
Gene Symbol	MANF
UniProt ID	MANF_HUMAN
Immunogen Sequence	Human
Specificity	MANF (Mesencephalic astrocyte-derived neurotrophic factor).

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---