

## Anti-Cerebral dopamine neurotrophic factor antibody (STJA0040777)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB
Host / Source	Mouse
Reactivity	Human

### PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	IgG
Dilution Range	WB 1:2000-1:4000
Formulation	PBS pH 7.4, with 0.1% sodium azide
Isotype	IgG1
Storage Instruction	Store at -20°C or -70°C upon receipt. After opening maintain undiluted in smaller aliquots at -20°C or -70°C for up to 6 months. Avoid multiple freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">441549</a>
Gene Symbol	<a href="#">CDNF</a>
UniProt ID	<a href="#">CDNF_HUMAN</a>
Immunogen Sequence	Human
Specificity	CDNF (conserved dopamine neurotrophic factor).

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	-------------------------------------------------------------------------------------------------------------------------