

## Anti-Pro-nerve growth factor antibody (Full-Length) [BS312] (STJA0040792)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | FC/ICC/WB          |
| Host / Source | Mouse              |
| Reactivity    | Human              |

### PRODUCT PROPERTIES

|                     |                                                                                                                                                                                                                                                            |
|---------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Clonality           | Monoclonal                                                                                                                                                                                                                                                 |
| Conjugation         | Unconjugated                                                                                                                                                                                                                                               |
| Purification        | IgG                                                                                                                                                                                                                                                        |
| Dilution Range      | WB 1-2 µg/mL<br>ICC 1-2 µg/mL                                                                                                                                                                                                                              |
| Formulation         | Lyophilized from a solution containing PBS buffer pH 7.4 with 3% trehalose, without preservatives.                                                                                                                                                         |
| Isotype             | IgG2b, lambda                                                                                                                                                                                                                                              |
| Storage Instruction | Spin vial briefly before opening. Reconstitute in 100 µL sterile-filtered, ultrapure water. Centrifuge to remove any insoluble material. Final buffer contains no preservatives. Store lyophilized antibody at 2-8°C. After reconstitution divide into ali |

### TARGET INFORMATION

|                    |                                            |
|--------------------|--------------------------------------------|
| Gene ID            | <a href="#">4803</a>                       |
| Gene Symbol        | <a href="#">NGF</a>                        |
| UniProt ID         | <a href="#">NGF_HUMAN</a>                  |
| Immunogen Region   | Full-Length                                |
| Immunogen Sequence | Human                                      |
| Specificity        | Human Species cross-reactivity not tested. |

### ADDITIONAL INFORMATION

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|